



Bionew™ Tissue DNA Purification Kit

Size: 100 preps

#Cat: BN100-012a

- I. Description
- II. Components and Storage Conditions
- III. Applicability
- IV. Materials Needed
- V. Tissue Pretreatment
- VI. Protocol
- VII. Troubleshooting
- VIII. Composition of Buffers and Solutions
- IX. Related Products

. Description

Nucleic acid isolation is the most technically demanding and laborious procedure performed in molecular cloning laboratory. At first purification procedure usually used phenol, chloroform, etc. These hazardous organic chemicals, need multiple centrifugations and will result in significant loss of DNA, may affect the enzyme digestion of DNA or introduce amplification inhibitors.

Using spin column to purify DNA from samples is conventional method now, which required DNA was adsorbed to the membrane by centrifuge or vacuum, which bring not only DNA but also cell debris and other precipitate to the membrane. The membrane can be clogged because of the precipitate, then result in inefficient recovery of samples. Especially for tissue sample, after grinding tissue sample and adding some lysis buffer to make a tissue lysate, the lysate is always containing some precipitate.

A recent developed DNA extraction method can directly use the lysate as the samples. It is the use of magnetic beads covered with nucleic acid-binding matrices. Its theory is nucleic acid (DNA and RNA) in the solution is specifically adsorbed to the surface of magnetic beads because the following action: (i) shielded intermolecular electrostatic forces, (ii) dehydration of the DNA/RNA and beads surfaces, and (iii) intermolecular hydrogen bond formation in the DNA/RNA-beads contact layer.

This method can farthest eliminate the inhibition in the DNA extraction course, because nucleic acid in the lysate is specifically adsorbed to the surfaces of magnetic beads but the precipitate in the lysate can't make any effect. Then after 2 times of wash steps for the surface of magnetic beads, the impurities are eliminated. Finally the pure DNA is eluted and can be used in PCR or enzyme digestion directly.

Bionew™ tissue DNA purification kit is the excellent kit which uses magnetic beads for DNA extraction. It can be widely used in DNA purification with sample from Research Institutions, Inspection and Quarantine, Life Science and Industrial Markets. The procedures of DNA extraction system contains: lysis and DNA release, DNA binding to surface of magnetic beads, 2 times of wash and DNA elution.

For more information, see the www.bionewtech.com or contact Bionewtech Technical Service: service@bionewtech.com



II. Components and Storage

Components	Size (100 preps)
Magnetic Beads Solution	800 μ l \times 1 tube
Proteinase K(Lyophilized)	2 tubes
Proteinase K Solution Buffer	1.4ml \times 2 tubes
Lysis Buffer	22 ml \times 1 bottle
Wash 1 Buffer	22 ml \times 1 bottle
Wash 2 Buffer	22 ml \times 1 bottle
Elution Buffer	11 ml \times 1 bottle

Upon receipt, store all components at room temperature until Expiry Date.

Expiry Date: From Production Date, for up to 12 months.

III. Applicability

The kit can be directly used for DNA purification from **tissue, bone, hair**, etc.

. Materials Needed

Reagents: 100% Alcohol , 100% Isopropyl Alcohol

Expendable: Centrifuge Tube of 2ml, 1.5ml ; Tip of 200 μ l, 1000 μ l ; All sterilized by 1.01E+5Pa, 121 for 20 min. Recommend to use Axygen expendable.

Instruments: Pipettor(10 μ l,200 μ l,1000 μ l); Vortex Agitator, Constant Temperature Equipment (water bath or metal bath)

. Tissue Pretreatment

- Weight 10mg tissue, adding liquid nitrogen to freeze tissue, and grind tissue to powder;
- After volatilization of liquid nitrogen, transfer tissue powder to 100 μ l [Lysis buffer] and mix equably, have a pretreatment of 56 ,30min;
- Centrifuge with 4000g, 2min, then transfer supernatant for next DNA purification;
- Pre-treat every tissue sample as above.

. Protocol

(1) Preparation Protocol

- 1.1 For the tube of [Proteinase K (Lyophilized)], take 1.4 ml [Proteinase K Solution Buffer] into the tube and solved the [Proteinase K (Lyophilized)] by up and down the tube. It is the final [Proteinase K]. Then make aliquots of [Proteinase K] and store at 2~8 until expiry date.
- 1.2 For the bottle of 22ml [Wash 1 Buffer], add 11ml Alcohol, 11ml Isopropyl alcohol to the bottle and mix it by the bottle up and down, mark it to confirm the addition of alcohol and isopropyl alcohol. It is the final [Wash 1 Buffer]. Screw the cap closely and store at 2~8 for up to 6 months.
- 1.3 For the bottle of 22ml [Wash 2 Buffer], add 22ml Isopropyl alcohol to the bottle and mix it by the bottle up and down, mark it to confirm the addition of isopropyl alcohol. It is the final [Wash 2 Buffer]. Screw the cap closely and store at 2~8 for up to 6 months.
- 1.4 Before use vortex [Magnetic Beads Solution] for 15s to mix equably. The dosage of beads is 7.5 μ l beads/sample, so for N samples take N \times 7.5 μ l beads into new centrifuge tube, then take N \times 200 μ l isopropyl alcohol into this tube, vortex the solution 15s to mix equably. It is the final [Magnetic Beads Solution].



- 1.5 If the [Lysis Buffer] forms a precipitate, warm the solution in 56 for 10mins to solve the precipitate completely.
- 1.6 For frozen samples take out to room temperature to unfreeze. Take out [Elution Buffer] in 56 to preheat it.

(2) Purification Protocol

- 2.1. For N samples, use N centrifuge tubes (2ml) and separately marked it; add 25μl step1.1 [Proteinase K] firstly.
- 2.2. Add 100μl tissue lysate from pretreated samples of **Part V**, using different 200μl tips
- 2.3. Add 100μl [Lysis Buffer] using 200μl tips, vortex for 5s, incubate the sample at 56 for 30min.
- 2.4. Add 200μl step1.4 [magnetic beads solution] using 200μl tips, vortex for 10s to mix it equably, incubate the sample at 56 for 3 minutes.
- 2.5. Place the N centrifuge tubes(2ml) containing sample mixture to the magnetic stand for 3 minutes, the magnetic beads are adsorbed to the wall of tube because of magnetic force, then discard all supernatant using different 1000μl tips.
- 2.6. Remove the tubes from the magnetic stand, add 400μl step1.2 [Wash 1 Buffer] using 1000μl tips, vortex 5s.
- 2.7. Place the tubes to the magnetic stand for 2 minutes, and discard all supernatant using 1000μl tips (permission of using the same 1000μl tip).
- 2.8. Remove the residual solution using 200μl tips.
- 2.9. Remove the tubes from the magnetic stand; add 400μl step1.3 [Wash 2 Buffer] into every tube with different 1000μl tips, mix every sample by slowly drawing and pumping by 1000μl pipettor. Note: Don't let tips touch to residual solution on the upside of tube.
- 2.10. Use N new centrifuge tubes(2ml), marked as step 2.1, transfer every mixture of step 2.9 into every new marked 2ml centrifuge tube
- 2.11. Place the tubes to the magnetic stand for 2 minutes, and discard all supernatant using 1000μl tips (permission of using the same 1000μl tip).
- 2.12. Open the lids of tubes and dry at room temperature for 3minutes, and then remove again the residual solution at the bottom of tubes with 200μl tips
- 2.13. Add 56 preheated [Elution Buffer] 100μl (Optional 50-100μl), mix every sample by slowly drawing and pumping by 200μl pipettor, incubate at 56 for 3 minutes.
- 2.14. Place the tubes to the magnetic stand for 2 minutes,
- 2.15. Transfer the supernatant to the new centrifuge tubes. It is the eluted genomic DNA, it can be used for next detection at once or stored at 4°C for 1 week or stored at -20 for 6 months.
- 2.16. The yield and purity can be determined by spectrophotometer. For double stand DNA:
dsDNA Concentration = $(OD_{260} - OD_{320}) \times 50 \mu\text{g/ml}$
dsDNA Purity = $(OD_{260} - OD_{320}) / (OD_{280} - OD_{320})$,



When using the frozen tissue as the sample, its concentration and purity is decreased comparing with fresh tissue, so the range of purity of pure dsDNA is between 1.7 and 2.0.

. Troubleshooting

For questions not addressed here, please contact Bionewtech.

Web site: www.bionewtech.cn

Email: service@bionewtech.cn

Symptoms	Possible Causes	Comments
Poor yield	Too little sample used , or the initial amount of DNA in the sample is poor.	Increase initial amount of sample.
	The tissue sample is frozen for a long time, and the DNA has been degraded	Use fresh samples or use more samples.
	Elution Buffer amount is not enough.	Be sure 100µl Elution is used to elute the DNA.
Beads in final elution	Beads is occasionally transferred by rapid pipetting	Vortex or mix solution, place in the magnetic stand and transfer eluate to new tube again.
Inconsistent yield	Inconsistent amounts of magnetic beads	Be sure to vortex the magnetic beads solution before use and vortex beads in the binding step.
	Inconsistent sample is used	Be sure to mix the sample equably before use.
Purity lower than 1.7	The tissue sample is frozen too long time.	If the purity is in the range of 1.6~1.7, the sample can be used for next detection. If the purity is lower than 1.6, don't recommend use such sample to do next detection.

. Composition of Buffers and Solutions

Elution Buffer: 10mM Tris. (PH7.2)

.Related Products

Bionew™ Whole Blood DNA Purification Kit	#Cat: BN100-011a
Bionew™ Tissue DNA Purification Kit	#Cat: BN100-012a
Bionew™ DNA Micro Purification Kit	#Cat: BN100-014a
Bionew™ Plasmid DNA Purification Kit (Medium Quantity)	#Cat: BN100-015a
Bionew™ Bacteria DNA Purification Kit	#Cat: BN100-016a
Bionew™ Viral DNA Purification Kit	#Cat: BN100-031a
Bionew™ Viral RNA Purification Kit	#Cat: BN100-032a
Bionew™ Viral DNA/RNA Purification Kit	#Cat: BN100-033a
Bionew™ Cells Total RNA Purification Kit	#Cat: BN100-021a
Bionew™ Tissue Total RNA Purification Kit	#Cat: BN100-022a
Bionew™ Whole Blood Total RNA Purification Kit	#Cat: BN100-023a
Bionew™ Lockable Magnetic Stand	#Cat: BN200-011a